

# Gator<sup>®</sup> High Sensitivity AAV9 Kit for AAV9 Titer

Gator<sup>®</sup> High Sensitivity (HS) AAV9 Kit is designed for AAV titer measurement in crude lysates and purified samples using Gator next-gen BLI systems. The titer for ultra crude samples can be accurately determined with the “dilute and dip” method using this kit. It is ideally suited for upstream and downstream monitoring, and final product QC. The dynamic range is comparable to ELISA and the titer values correlate well. In addition, the kit offers significant benefits such as automation and faster analysis.

## PRODUCT INFORMATION

**Part Number**  
350005

### Includes

- HS AAV9 probes (96 probes/tray)
- HS AAV9 Detection Solution
- AAV Amplification Solution
- AAV Substrate Solution
- AAV Substrate Diluent
- Q Buffer

## PERFORMANCE SUMMARY

**Dynamic Range**  
1 x 10<sup>7</sup> - 1 x 10<sup>9</sup> vp/mL


**Assay Time**  
8 samples in 35 minutes  
24 samples in 99 minutes

**Crude Sample Tolerant**  
Yes

## KEY BENEFITS

- Sensitivity comparable to ELISA
- Crude sample tolerant
- Dilute and dip
- Fully automated including washes

## DYNAMIC RANGE



AAV9 standard curve generated using the HS AAV9 kit. The concentration is based on end point nanometer shift.

## ACCURACY AND PRECISION

Conc (vp/mL)	nm Shift	Calc Conc (vp/mL)	% Recovery	% CV (n=3)
1.00E+09	12.7	1.02E+09	102	1.7
3.33E+08	4.31	3.30E+08	99	4.4
1.11E+08	1.48	1.09E+08	98	2.3
3.70E+07	0.49	3.48E+07	94	2.6
1.23E+07	0.18	1.24E+07	101	4.4

Percent recovery for HS AAV9 kit is above 90% in Q Buffer. The kit shows very good precision with a CV below 10%.



**CRUDE SAMPLES USING DILUTE AND DIP**

<b>Conc (vp/mL)</b>	<b>Calc Conc (vp/mL)</b>	<b>% Recovery</b>	<b>% CV (n=3)</b>
1.00E+09	1.09E+09	109	4.3
1.11E+08	1.01E+08	93	5.7
1.23E+07	1.26E+07	102	9.7

The HS AAV9 Kit is compatible with crude matrices such as cell lysates and media. The table shows recovery and CV for AAV9 in HEK 293-T cell lysate diluted 1:10 in Q Buffer. AAV9 capsid was purchased from [www.virovek.com](http://www.virovek.com).

