

## Bisulfite Sequencing with the DS3000 Compact CE Sequencer

### Abstract


Of the four bases that make up genomic DNA in higher eukaryotes, the fifth carbon of cytosine can be methylated by a methyltransferase. Cytosine methylation plays an important role in various events that occur on genomic DNA, including the regulation of gene expression, and it is known as a representative epigenetic marker. Bisulfite sequencing is a commonly accepted method of analyzing methylated cytosine on a single-nucleotide basis. In bisulfite sequencing, extracted genomic DNA is treated with bisulfite (hydrogen sulfite) for a certain period. During this time, cytosine is deaminated and converted to uracil, while methylated cytosine remains as it is. Therefore, by sequencing the bisulfite-treated DNA, it is possible to distinguish between cytosine and methylated cytosine. Here, we introduce a practical example of bisulfite sequencing using the DS3000 Compact CE Sequencer (DS3000). We used genomic DNA extracted from the HCT116 cell line, which is artificially methylated by CpG methylase. The bisulfite-treated sample was subjected to electrophoresis on DS3000, and the sequences before and after bisulfite treatment were compared. As a result, 19 out of 20 methylated cytosines were detected.



**DS3000**  
Compact CE Sequencer


### Results

First, HCT116 genomic DNA was treated with bisulfite, then purified and quantified. Unmethylated cytosines were converted to uracil in this step. Next, the promoter site of the human *MLH1* gene was amplified by polymerase chain reaction (PCR). PCR amplification converts uracil to thymine, while methylated cytosine remains as cytosine. *MLH1* is involved in mismatch repair and its promoter site is strongly methylated in HCT116 cells. Next, the PCR products were treated with the BigDye<sup>®</sup> Terminator v3.1 Cycle Sequencing Kit (Applied Biosystems<sup>®</sup>). Finally, each reaction product was analyzed by DS3000. Representative results are shown in Figure 1. Sharp and clear peaks were observed throughout the entire window of the electropherogram. A mixed base (Y) of cytosine (C) and thymine (T) can be seen at the 7<sup>th</sup> base of the electropherogram. This is presumably due to the residual cytosine in some DNA molecules. The positions of the detected methylated cytosines are summarized in Figure 2. In the *MLH1* promoter, 20 CpG sequences that serve as substrates for CpG methylase are known to be present. Among these sequences, peaks derived from methylated cytosine were identified in 19 sequences, excluding the mixed base (Y) mentioned above.



**Figure 1: Representative electropherogram of bisulfite-treated DNA**

Results obtained from the forward primer are shown. Please refer to the main text for details on the mixed base (Y).



**Figure 2: Locations of the methylated cytosines**

The observed methylation pattern was determined using the contig sequence from the forward and reverse sequences.

The first 40 bases of each sequence were excluded. Before bisulfite: before bisulfite treatment. After bisulfite: after bisulfite treatment.

**Conclusion**

Here, we described a practical example of bisulfite sequencing using DS3000. By targeting the human *MLH1* gene promoter as a model system, 19 out of 20 methylation sites were detected. The remaining undetected site was thought to be due to residual cytosine molecules from incomplete bisulfite treatment. DS3000 is equipped with two electrophoresis protocols: fast sequencing and standard sequencing. The two protocols provided results of equivalent quality.

**Troubleshooting**

Bisulfite-treated DNA molecules become single-stranded and are unstable. We recommend handling them with extreme care. Keep the sample on ice whenever possible, and shorten the handling time.

## Materials and Methods

First, 500 ng of EpiScope<sup>®</sup> Methylated HCT116 g DNA (TaKaRa<sup>®</sup>) was treated with bisulfite using the EZ DNA Methylation-Startup<sup>™</sup> Kit (Zymo Research). After the treatment, the DNA was purified, eluted in 10  $\mu$ L of elution buffer, and 1  $\mu$ L of the solution was subjected to PCR. PCR was performed with the conditions summarized in Tables 1 to 3 using control primer included in the kit. After purification with DNA Clean & Concentrator<sup>®</sup>-5 (Zymo Research), the cycle sequencing reaction was performed using 200 pg of PCR product and the BigDye<sup>®</sup> Terminator v3.1 Cycle Sequencing Kit (Applied Biosystems<sup>®</sup>). The resulting product was purified by ethanol precipitation, and subjected to capillary electrophoresis using DS3000 with both the fast and standard sequencing protocols (Secondary Peak Height Threshold was set to 25%). The electropherogram was visualized by SangerseqR (1) with modifications.

**Table 1: PCR cycle**

| Temperature(°C) | Time (s) | Cycle |
|-----------------|----------|-------|
| 94              | 600      | ×1    |
| 94              | 30       | ×35   |
| 56              | 30       |       |
| 72              | 60       |       |
| 72              | 420      | ×1    |
| 4               | hold     |       |

**Table 2: PCR mixture**

| Reagent                           | Volume ( $\mu$ L) |
|-----------------------------------|-------------------|
| 2× ZymoTaq <sup>™</sup> PreMix    | 12.5              |
| Forward primer (10 pmol/ $\mu$ L) | 2                 |
| Reverse primer (10 pmol/ $\mu$ L) | 2                 |
| DNA (~10 ng/ $\mu$ L)             | 1                 |
| PCR-grade water                   | 7.5               |

**Table 3: PCR primer sequences**

| Primer  | Sequence                        |
|---------|---------------------------------|
| Forward | GGAGTGAAGGAGGTTACGGGTAAGT       |
| Reverse | AAAAACGATAAAACCCTATACCTAATCTATC |

### Reference

1, Hill, J., T., et al. Poly peak parser: Method and software for identification of unknown indels using sanger sequencing of polymerase chain reaction products. *Developmental Dynamics*. **243**, 1632-1636 (2014)

## Main specifications for the DS3000

### Main unit specifications

| Item                                  | Details  |
|---------------------------------------|--|
| Number of capillaries                 | 4  |
| Capillary length                      | 36 cm  |
| Sample format                         | 8-tube strip × 4   |
| Device control                        | Touch panel PC   |
| Number of Dyes                        | 6  |
| Application                           | Sequencing analysis / Fragment analysis  |
| Size                                  | 400 (W) × 600 (D) × 600 (H) mm   |
| Weight                                | 45 kg  |
| Performance guarantee temperature     | 15 – 30°C  |
| Performance guarantee humidity        | 20 – 80% RH (no condensation)  |
| Power input                           | 100 – 240 ±10% VAC, 50/60 Hz   |
| Rated power                           | 260 VA   |
| Supported secondary analysis software | <ul style="list-style-type: none"> <li>• Mutation Surveyor (SoftGenetics, LLC, sold separately)</li> <li>• GeneMarker (SoftGenetics, LLC, sold separately)</li> <li>• GeneMarker HID (SoftGenetics, LLC, sold separately)</li> </ul> |

### Run module specifications

| Run Module                       | Application             | Polymer type | Contiguous Read Length* <sup>1</sup> (bp, QV20 CRL) | Average run time (minutes) |
|----------------------------------|-------------------------|--------------|---|----------------------------|
| Fast_Sequence36_Polymer7         | Sequencing analysis     | Polymer7     | ≥600  | ≤32                        |
| Standard_Sequence36_Polymer7     | Sequencing analysis     | Polymer7     | ≥700  | ≤60                        |
| BDx_Fast_Sequence36_Polymer7     | BDx sequencing analysis | Polymer7     | ≥600  | ≤32                        |
| BDx_Standard_Sequence36_Polymer7 | BDx sequencing analysis | Polymer7     | ≥700  | ≤60                        |

| Run Module                   | Application       | Polymer type | Average run time (minutes) | Sizing precision* <sup>2</sup> (bp, 50 – 400 bp) |
|------------------------------|-------------------|--------------|----------------------------|--|
| Fragment_Analysis36_Polymer7 | Fragment analysis | Polymer7     | ≤35                        | NA   |
| Fragment_Analysis36_Polymer4 | Fragment analysis | Polymer4     | ≤44                        | <0.16  |

\*1 Contiguous Read Length (bp, QV20 CRL) is measured with BigDye® Terminator v3.1. Sequencing Standard Kit (Thermo Fisher Scientific, sold separately)

\*2 Sizing precision (bp, 50-400 bp) is measured with PowerPlex™ ESI17 Fast Allelic Ladder and WEN ILS 500 ESS. (Promega®, sold separately)

• BigDye is a registered trademark of Thermo Fisher Scientific Inc.

• Promega, PowerPlex is a registered trademark of Promega Corporation.

### Consumables specifications

| Product name                          | Part number | Details  | Remarks                        |
|---------------------------------------|-------------|--|--------------------------------|
| Capillary Cartridge 36 cm             | 613-0330    | 1 pcs  | Storage temperature: 15 – 30°C |
| Buffer                                | 613-0252    | Anode Buffer × 2 cartridges<br>Cathode Buffer × 2 cartridges | Storage temperature: 2 – 10°C  |
| Polymer7                              | 613-0251    | 4 cartridges   | Storage temperature: 2 – 10°C  |
| Polymer4                              | 613-0250    | 4 cartridges   | Storage temperature: 2 – 10°C  |
| Septa for Cathode Buffer Cartridge    | 613-7231    | 10 pcs   |                                |
| Retainer for Cathode Buffer Cartridge | 613-7233    | 4 pcs  |                                |
| Septa for 8 well tubes                | 613-7230    | 24 pcs   |                                |
| Base and Retainer for 8 well tubes    | 613-7232    | 4 pcs  |                                |
| Anode Electrode Assembly              | 613-7263    | 1 pcs  |                                |

Specifications in this catalog are subject to change with or without notice, as Hitachi High-Tech Corporation continues to develop the latest technologies and product for its customers.

CAUTION:For correct operation, follow the instruction manual when using the instrument.

NOTICE:The system is For Research Use Only, and is not intended for any animal or human therapeutic or diagnostic use. Hitachi High-Tech does not guarantee the performance on the document with every possible sample under every possible condition. Copyright (C) Hitachi High-Tech Corporation 2022 All rights reserved.

 Hitachi High-Tech Corporation

Tokyo, Japan

[www.hitachi-hightech.com/global/science/](http://www.hitachi-hightech.com/global/science/)

Toranomon Hills Business Tower, 1-17-1 Toranomon, Minato-ku, Tokyo 105-6409, Japan  
customercenter.ev@hitachi-hightech.com

